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GB TABLE OF CONTENTS

1	For your safety and that of your patients 35		5	Installation and operation..... 42
1.1	Symbols and information on the appliance35		5.1	Connection to a compressed gas cylinder 43
1.2	General safety information 35		5.2	Connection/Assembly of downstream equipment46
1.3	Safety information for working with medical gases37		5.3	Commissioning47
1.4	Rules and guidelines 38		5.4	At the end of operation..... 49
2	Product description 38		5.5	Changing the compressed gas cylinder50
2.1	Intended purpose38		5.6	Decommissioning and storage..... 50
2.2	Special regulations for pressure regulators with pressure outlet39		6	Maintenance and care .. 51
2.3	Special regulations for pressure regulators with flow outlet39		6.1	Cleaning and disinfecting 51
2.4	Special regulations for pressure regulators with female thread connection..... 39		6.2	Performance test52
2.5	Special regulations for pressure regulators for test gas mixtures.....39		6.3	Maintenance and service life54
2.6	Special regulations for stainless steel pressure regulators..... 40		6.4	Warranty, customer service and dispatch 55
2.7	Special regulations for AMAG pressure regulators.....40		6.5	Troubleshooting56
2.8	Special regulations for pressure regulators40		6.6	Disposal56
2.9	Set-up and mode of operation40		7	Combination with downstream equipment .57
3	About the instructions for use..... 41		8	Appendix 57
3.1	Scope.....41		8.1	Pressure regulator: Technical specifications.....57
3.2	Target audience41		8.2	Pressure relief valve: Technical specifications.....59
3.3	Improvements42		8.3	Permissible ambient conditions59
4	Transport and packaging 42		8.4	Flow curves (for AIR)60
			8.5	Classification61
			8.6	Classification / UMDNS code.....61
			8.7	Ordering information61
			9	Special information61
			9.1	Information on the electronic high-pressure sensor61

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




3-Leitersystem / 3-wire system

⊖ → : Hilfsenergie/Power Supply

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
1 FOR YOUR SAFETY AND THAT OF YOUR PATIENTS

1.1 Symbols and information on the appliance

The following section offers an explanation of the meaning of the symbols that are both attached to the appliance (e.g. type plate) and contained in or on the accompanying information and packaging.

GB

	<p>"Caution"</p> <p>This symbol indicates instructions which, in the event of non-compliance, could result in injury and/or property damage.</p>
	<p>"Beware"</p> <p>This symbol indicates instructions which, in the event of non-compliance, could result in property damage.</p>
	<p>"Free from oil and grease"</p> <p>This symbol denotes equipment that uses medical gases and must therefore be kept free from oil and grease.</p>
	<p>"Manufacturer"</p> <p>This symbol denotes the name and, possibly, address of the manufacturer.</p>
	<p>"Date of manufacture"</p> <p>This symbol denotes the date of manufacture.</p> <p>XXXX - XX → Month of manufacture ↓ Year of manufacture</p>
	<p>"Temperature restriction"</p> <p>This symbol indicates a temperature range, e.g. in connection with storage or operating conditions.</p>
	<p>"Relative humidity limits"</p> <p>This symbol indicates the range of relative humidity that is permitted, e.g. for storage or operating conditions.</p>
	<p>"Ambient pressure limits"</p> <p>This symbol indicates the range of ambient pressure that is permitted, e.g. for storage or operating conditions.</p>

1.2 General safety information

Please read and observe the following information before using or performing maintenance on the TESCOM EUROPE pressure regulator.

GB

- The appliance may only be used for its intended purpose as described in these instructions for use (see Section "Intended purpose").
- These instructions for use are integral to the pressure regulator and must be available to the user at all times.
- You must understand how the pressure regulator works before commencing operation with it.
- Never use an appliance if it is leaking or faulty!
- Operation of the pressure regulator enriches and, where appropriate, controls the corresponding gas in the patient environment. The pressure regulator may therefore only be used in accordance with a doctor's instructions.
- Medical gases are a medicinal product and may only be administered under the supervision of qualified medical staff who are familiar with the precautions that need to be observed for the respective treatment.
- The safety data sheets and information provided by the gas manufacturers must be observed.
- Before commissioning the pressure regulator, make absolutely sure that the gas type, the gas volume flow range and the connection are correct and that any upstream or downstream equipment has been fitted securely.
- Not suitable as a drive for surgical equipment.
- Only accessories bearing the CE mark that are compatible with the pressure regulator both in terms of intended use and connection may be used (see Section "Combination with downstream equipment").
- To avoid misunderstandings, never leave the appliance unattended or connected to downstream equipment when it is not in use.
- Never try to open the appliance since this could result in the release of spring-loaded parts which might cause injury.
- To avoid the risk of infection or bacterial contamination, please read and observe the information in the Section "Maintenance and care".
- The pressure regulator may not be immersed in liquid.
- The pressure regulator may not be sterilised.
- Foreign materials or fluids must not get inside the appliance. This applies to dust and disinfectants in particular. Use protective caps / plugs, either those supplied or other suitable ones in clean condition.
- You **MUST** wash your hands before working on the gas supply.
- Ensure that the appliance cannot come into contact with oil, grease or any other lubricants during operation or storage. Oil and grease, as found for example in hand creams, can react violently with pressurised gases (especially oxygen and nitrous oxide). Risk of explosion!
- Smoking and naked flames are strictly forbidden in the vicinity of valves and fittings carrying gas.
- Always proceed with caution when opening the cylinder valve.
- The pressure regulator connections must match those on the gas supply or downstream equipment. Adapters must not be used.
- The compressed gas cylinder must be secured at all times so that it cannot tip over.

- Make sure that the temperature of the gas administered corresponds to the patient's body temperature. Risk of injury to the patient.
 - Do not replace or remove the filter.
 - Do not block or cover the pressure relief valve.
 - Do not perform any technical or design modifications.
 - Comply with the storage and operation conditions.
 - Do not use the appliance for direct supply to the patient.

- ➔ • Do not apply any external loads or forces; support or secure any connected appliances.

- Do not hang objects on it.

- The pressure regulator must not be used as a handle on mobile appliances or during transport with the gas cylinder connected.

- Protect against impacts, jolts, falls and vibrations.

- Special safety equipment must be used on appliances connected to the pressure regulator to ensure that no dangerous pressure can accumulate inside them. The pressure relief valve on the pressure regulator is not suitable for protecting these appliances.

1.3 Safety information for working with medical gases

Observe the gas manufacturers' relevant safety information and instructions for use.

Oxygen (O₂) increases combustibility and flammability as well as encouraging flammable materials to burn. Greases and oils (such as those found in hand creams or lubricants) can ignite on contact with pressurised oxygen. Good ventilation is absolutely essential. High levels of concentration can be damaging to health. Smoking and naked flames are forbidden.

Nitrous oxide (N₂O) has a strong anaesthetic effect. It increases combustibility and flammability as well as encouraging flammable materials to burn. Greases and oils (such as those found in hand creams or lubricants) can ignite on contact with pressurised nitrous oxide. High levels of concentration could cause anoxia and apnoea. Good ventilation is absolutely essential. Smoking and naked flames are forbidden.

Carbon dioxide (CO₂) has an intoxicating effect at low concentrations, and a debilitating effect (apnoea) at higher concentrations. Good ventilation is absolutely essential.

Xenon (Xe) has a strong anaesthetic effect. High levels of concentration could cause anoxia and apnoea. Good ventilation is absolutely essential.

Test gas mixtures consist of a mixture of gases, some of which are dangerous. In certain compositions, and depending on the pressure and temperature of the test gas mixtures and the valve materials, dangerous reaction and/or decomposition products may occur. Prior to using the test gas mixtures, it is absolutely essential that you read the gas manufacturers' relevant safety information and instructions for use and that you ensure that the valves and fittings you are using are compatible for the relevant gases.

1.4 Rules and guidelines

Compliance with national laws, regulations and guidelines on the safe handling and operation of medical equipment in the context of 2007/47/EC and valves and fittings that carry gas, on handling gas under pressure, especially oxygen (O₂) and nitrous oxide (N₂O), with the information provided by the suppliers and manufacturers of the compressed gases and with the instructions for using any connected equipment is mandatory.

2 PRODUCT DESCRIPTION

2.1 Intended purpose

The pressure regulator is used for a variety of medical equipment in the field of medicine and therapy to regulate the pressure at which medical gases and/or gas mixtures are administered.

- For connection by a trained user to a compressed gas cylinder.
- Only for use with the medical gas and/or gas mixture indicated on the pressure regulator.
- Medical gases and gas mixtures are a medicinal product and may only be administered under the supervision of qualified medical staff who are familiar with the precautions that need to be observed for the respective treatment. The safety data sheets and information provided by the gas manufacturers must be observed.
- Must not be used as a drive for surgical equipment.
- Not suitable for operation in the range of strong magnetic fields, e.g. those generated by magnetic resonance imaging, unless otherwise indicated.
- For stationary and mobile use, land and air transportation.
- For supplying suitable medical equipment in accordance with the technical specifications indicated on the type plate and in the instructions for use.
- The TESCOM EUROPE pressure regulator is not fitted with an output pressure gauge. The connected appliance must be used if you wish to monitor the output pressure during operation.

Certain versions of the pressure regulator are subject to additional, and possibly complementary, special regulations:

2.2 Special regulations for pressure regulators with pressure outlet

- Reduction of 30 MPa (300 bar) inlet pressure maximum to a nominal outlet pressure of 400 kPa (4 bar).
- Not suitable for use as a flow control element.

2.3 Special regulations for pressure regulators with flow outlet

- Reduction of 30 MPa (300 bar) inlet pressure maximum to a maximum outlet pressure of 550 kPa (5.5 bar).
- For setting the gas flow to the patient's needs.
- For inhalation or insufflation to the nose and/or throat area for patients with stable spontaneous respiration. The gas pressure is not sufficient to overcome the resistance of the respiratory tract.

- For flushing cavities with gas using additional suitable measures to constantly monitor and control the output pressure on the appliance. The pressure in the cavity can rise to up to approx. 5.5 bar (> 4,100 mm Hg).

2.4 Special regulations for pressure regulators with female thread connection

- Reduction of 30 MPa (300 bar) inlet pressure maximum to a nominal outlet pressure p_2 at normal flow Q in accordance with the type plate on the pressure regulator.
- This pressure regulator may only be operated as a fixed component of medical appliances or must be permanently (i.e. not detachable without tools) connected to medical equipment.
- Not suitable for use as a flow control element.

2.5 Special regulations for pressure regulators for test gas mixtures

- Reduction of 20 MPa (200 bar) inlet pressure maximum to a nominal outlet pressure p_2 at normal flow Q in accordance with the type plate on the pressure regulator.
- Only suitable for gas mixtures with a lower oxidising effect than air.
- Depending on the respective composition of the test gas mixture, the pressure, the temperature of the gas and the time spent by the test gas mixture in the pressure regulator and in any connected equipment, dangerous, e.g. very toxic and/or flammable, products may occur. Prior to use, the user must always carefully check to ensure that the relevant test gas mixture is compatible with the pressure regulator and connected equipment. Please do not hesitate to contact TESCOM EUROPE if you require any information about the materials and surface coatings of the pressure regulator. You will find the contact address at the end of these instructions for use.
- Only for stationary use in hospitals and clinics.

2.6 Special regulations for stainless steel pressure regulators

GB

- Reduction of 20 MPa (200 bar) inlet pressure maximum to a nominal outlet pressure p_2 at normal flow Q in accordance with the type plate on the pressure regulator.
- Only suitable for gas mixtures with a lower oxidising effect than air.
- Only suitable for operation at gas and ambient temperatures of between 0 °C and 60 °C.
- Only for stationary use in hospitals and clinics.

2.7 Special regulations for AMAG pressure regulators

Reduction of 20 MPa (200 bar) inlet pressure maximum to a nominal outlet pressure of 400 kPa (4 bar).

The WegaMed AMAG pressure regulator is suitable for use in the range of magnetic fields up to 3 Tesla (e.g. magnetic resonance imaging, etc.).

2.8 Special regulations for pressure regulators

2.9 Set-up and mode of operation

The TESCOM EUROPE pressure regulator is a single-stage, direct-acting, spring-loaded pressure regulator without primary pressure compensation. A diaphragm made of fibre-reinforced elastomer acts as the seal to the atmosphere.

The pressure regulator is fitted with an internal pressure relief valve. In the event of malfunction, medical gas can flow in the area surrounding the pressure regulator. Observe the safety information for working with medical gases.

Some versions are equipped with a locking flow outlet for adjusting the volume flow of the gas. A hand wheel is used to position an orifice plate containing different orifices to release the set flow. The possible gradations are specified by the appliance. Settings between these gradations are not possible.

Some versions are equipped with a continually adjustable flow outlet for measuring and adjusting the volume flow of the gas. The gas flows upwards against the weight of the float. The float adjusts to a certain height, depending on the volume flow. The volume flow

is adjusted by a needle valve downstream of the flow gauge.

Some versions of the pressure regulator are equipped with an alarm feature on the input connection. An alarm sounds if the inlet pressure drops below a specific value (see marking on the alarm feature). Medical gas flows into the area surrounding the pressure regulator during this process. Observe the safety information for working with medical gases.

The materials in the pressure regulator are compatible with all commonly used medical gases and/or mixtures of the same which contain oxygen, nitrous oxide, carbon dioxide and xenon and/or helium. If other gas components are used, their compatibility with the materials must be checked in each individual case.

3 ABOUT THE INSTRUCTIONS FOR USE

3.1 Scope

These instructions for use are valid for the TESCOM EUROPE WegaMed with the part number indicated on the last page of these instructions for use.

Certain versions may be subject to further instructions for use, in addition to these.

3.2 Target audience

- These instructions for use are intended for anyone responsible for the installation, operation or monitoring of the above valves and fittings.
- Maintenance may only be performed by TESCOM EUROPE or companies authorised by TESCOM EUROPE.

- The operator has a duty to ensure that only personnel who are familiar with the general practice and safety measures associated with working with gases and gas valves and fittings, as well as with the fundamental regulations with regard to occupational safety and accident prevention, are permitted to work on the device. These personnel must also have been taught how to handle this pressure regulator, and have read and understood the instructions for use.

- The design of the appliance, the accompanying documentation and the marking are based on the assumption that the appliance will only be installed, operated and monitored by trained and qualified personnel and that the user is familiar with certain general characteristics of the appliance. As such, the instructions, warnings and safety information focus mainly on features that are particular to the TESCOM EUROPE pressure regulator.

- The accompanying documentation does not address various general hazards and risks that should be obvious to qualified medical personnel. These include the improper use of the equipment and the possible side effects on patients in an unusual state of health.

- The user of the pressure regulator is responsible for taking all of the necessary safety precautions to ensure the safe operation of the equipment. The user and operator are responsible for deciding which precautions to take.
- These instructions for use form an integral part of the pressure regulator. Operators must be given access to these instructions. Additional copies can be supplied on request.

3.3 Improvements

If you have any suggestions for improving the product or the instructions for use, we would be delighted to hear from you. Both our products and the instructions for use are subject to a process of continuous development and improvement. The address and telephone number for TESCOM EUROPE GMBH & CO. KG appear on the last page of these instructions for use.

4 TRANSPORT AND PACKAGING


The appliance may only be transported (e.g. delivery to the manufacturer) and stored in suitable, sturdy packaging (e.g. the original packaging). During transport and storage the inlet and outlet connections must be covered by suitable and clean protective caps/plugs.

5 INSTALLATION AND OPERATION

- The user and operator are responsible for using and operating the appliance and accessories.
- You must wash your hands thoroughly before commencing any work involving the gas supply. Hydrocarbon compounds (e.g. oils, greases, rubbing alcohols, hand creams or remains of adhesive plasters) can cause explosive reactions if they come into contact with highly compressed gases.
 -  Medical gases and gas mixtures are a medicinal product and may only be administered under the supervision of qualified medical staff who are familiar with the precautions that need to be observed for the respective treatment. The safety data sheets and information provided by the gas manufacturers must be observed.

5.1 Connection to a compressed gas cylinder

The use of dust-free, dry medical gases is essential if the TESCOM EUROPE pressure regulator is to function as intended.

- | | | | |
|-----|---------------------------------|-----|---|
| 1 | Basic appliance | 2.5 | Gasket for pin index |
| 2 | Inlet connection | 2.6 | Gasket for threaded connection
(profile gasket or flat sealing disc) |
| 2.1 | Handwheel connection | 3 | Outlet connection
(specific to the type of gas) |
| 2.2 | Pin index | 4 | Gauge |
| 2.3 | Threaded connection | | |
| 2.4 | Gasket for handwheel connection | | |

Medical gases and gas mixtures are a medicinal product and may only be administered under the supervision of qualified medical staff who are familiar with the precautions that need to be observed for the respective treatment. The safety data sheets and information provided by the gas manufacturers must be observed.

Make sure that there is no evidence of damage or soiling on the connections on the compressed gas cylinder and the pressure regulator.

Damaged gaskets must be replaced. New gaskets can be purchased from TESCOM EUROPE or an authorised specialist dealer.

1. You must wash your hands thoroughly before commencing any work involving the gas supply. Hydrocarbon compounds (e.g. oils, greases, rubbing alcohols, hand creams or adhesive plaster) can cause explosive reactions if they come into contact with highly compressed gases.

- Using the technical specifications on the type plate and in these instructions for use, check whether the pressure regulator in your possession is technically suitable for the intended purpose (type of gas, pressure, flow). The maximum permissible inlet pressure of the pressure regulator must be greater than or equal to the filling pressure of the compressed gas cylinder. Read the information provided by the supplier of the compressed gases carefully.

Only use amagnetic compressed gas cylinders in the range of strong magnetic fields. Check that they are suitable prior to use.

- Remove any protective caps, seals, fasteners, etc. from the outlet connection of the cylinder valve.
- Before connecting the TESCOM EUROPE pressure regulator, slowly open the valve on the compressed gas cylinder briefly in a well ventilated room or outdoors (this is to blow out any contamination).

Do not aim the gas jet at anyone.
Please also observe the Section "Safety information for working with medical gases".

- Remove the protective caps covering the connections on the pressure regulator and keep them in a safe place.
- Screw the TESCOM EUROPE pressure regulator onto the compressed gas cylinder. The connections must match exactly – **do not use adapters!**

All connections must be clean and free from oil and grease. Do not use lubricants.

Connection using the handwheel

- The pressure regulator can be fitted with a handwheel connection to make it easier to connect it to the cylinder valve. This connection is denoted by the knurled connection screw or nut and the elastomer gasket inside one of the grooves on the coupling connection.

When making the connection, it is absolutely essential to check the gasket for evidence of damage and soiling. If the gasket is damaged or soiled, the pressure regulator must not be used and the gasket has to be replaced (see Section "Replacing gaskets").

Connection with pin index

- The pressure regulator can also be fitted with a pin index. Please note that only the original TESCOM EUROPE gasket included in the scope of supply may be used. When making the connection, it is absolutely essential to check the gasket for evidence of damage and soiling. If the gasket is damaged or soiled, the pressure regulator must not be used and the gasket has to be replaced (see Section "Replacing gaskets").

- The pins sticking out of the pin index must fit into the bore on the cylinder valve. The arrangement of these pins will vary depending upon the type of gas.

Connection by means of a threaded connection


- The pressure regulator may be fitted with a threaded connection. In this case, the connection screw takes the form of a hexagon nut with a right-handed or left-handed thread. Nuts with left-handed threads are denoted by a groove all the way around them.
- Threaded connections with right-handed threads are made by rotating the connection screw clockwise onto the cylinder valve.
- Threaded connections with left-handed threads are made by rotating the connection screw anticlockwise onto the cylinder valve.
- Threaded connections may only be tightened using an appropriate tool.
- When making threaded connections, never use the pressure regulator, manometer or outlet connection as a means of securing them.

The gasket may take the form of a profile gasket or flat sealing disc. When making the connection, it is absolutely essential to check the gasket for evidence of damage and soiling. If the gasket is damaged or soiled, the pressure regulator must not be used and the gasket has to be replaced (see Section "Replacing gaskets").

- The threaded connection must be made on the cylinder valve so that it is hand-tight. When you are doing this, align the pressure regulator as appropriate.
- When tightening the hexagon nut, the hexagon or flat face on the connection piece of the pressure regulator must be used to secure it.

Fig. 2

7. Check that the pressure regulator and the compressed gas cylinder are securely connected before slowly opening the cylinder valve. If the connection is not tight, close the cylinder valve immediately.

GB

Pressure surges caused by opening the valve too quickly can cause explosions!

8. A performance test must be run on the pressure regulator after every cylinder change. Should faults be detected during the performance test, the pressure regulator must be taken out of service. Contact TESCO EUROPE or an authorised specialist dealer to have the fault rectified.

A complete performance test includes:

- A visual inspection for mechanical damage
- Checking the system for leaks
- Checking the pressure relief valve for leaks

The user and operator are responsible for using and operating the pressure regulator and the equipment connected to it, together with any accessories.

5.2 Connection/Assembly of downstream equipment

A variety of equipment can be connected downstream of the TESCO EUROPE pressure regulator.

Only suitable equipment may be connected to the pressure regulator. The Section "Combination with downstream equipment" describes which accessories are permissible.

When connecting equipment, please note the following:

- The pressure regulator, the equipment to be connected and the joining elements must be performance tested and the compatibility of the materials with the gas and/or gas mixture verified.
- All downstream appliances, hoses and joining elements must be marked with the same type of gas.
- Connections and gaskets must be clean and free from oil and grease.
- There must be no evidence of damage on the gasket and sealing surfaces.
- Correct and proper procedures must be followed when connecting the TESCO EUROPE pressure regulator to downstream equipment. Please also read the instructions issued by the manufacturers of the downstream equipment you are using.

- Only use amagnetic accessories and appliances in the range of strong magnetic fields. Check that they are suitable prior to use.

- The connections must match - do not use adapters!
 - The opening of the pressure relief valve must not be covered or blocked.
 - Once the relevant connections have been made, a check must be carried out to ensure that all equipment is in working order. Please also follow the instructions for use issued by the manufacturers of the downstream equipment you are using.
 - The connection must be gas tight and must be checked once assembly is complete.
- ➔
- Do not apply any external loads or forces to the connections; support and/or secure any connected appliances and accessories.
 - Protect against impacts, jolts, falls and vibrations.

5.3 Commissioning

Prior to using the pressure regulator, briefly open and close the valve on the compressed gas cylinder twice in well ventilated rooms or outdoors in order to flush the pressure regulator.

Do not aim the gas jet at anyone.
Please also observe the Section "Safety information for working with medical gases".

Connect and start up downstream equipment. Please also follow the instructions for use issued by the manufacturers of the downstream equipment you are using.

Medical gases and gas mixtures are a medicinal product and may only be administered under the supervision of qualified medical staff who are familiar with the precautions that need to be observed for the respective treatment. The safety data sheets and information provided by the gas manufacturers must be observed.

Open the valve on the compressed gas cylinder slowly. Open the handwheel on the cylinder valve as far as it will go and then turn it back half a revolution. The cylinder pressure will now be displayed on the inlet pressure gauge.

Pressure surges caused by opening the valve too quickly can cause explosions!

Take suitable precautions to ensure that gas from the patient is not able to get inside the pressure regulator (backflow). Should this happen, the appliance must be sent back to the manufacturer immediately for cleaning.

The temperature of the gas and the pressure regulator drops sharply during operation. This can produce hypothermia in the patient.

GB

The pressure in the connected compressed gas cylinder and the patient must be monitored constantly during operation.

If no gas is being drawn: close the valve on the compressed gas cylinder. The pressure regulator is not a shutoff valve.

5.3.1 Flow outlet

The following additional instructions apply when commissioning versions with flow outlet.

- Turn the handwheel on the flow outlet clockwise to "0".
- Immediately before starting to treat the patient, all gas supply equipment (including downstream equipment) must be flushed at max. gas volume flow. This will ensure that any gas in the system is replaced. The flushing duration is determined by the downstream equipment.
- Turn the handwheel to select the required gas volume flow. The handwheel is turned anticlockwise to increase the flow.

5.3.1.1 Locking adjustable flow outlet

A stop prevents the gas volume flow being switched directly between zero and maximum and vice versa.

The handwheel must engage securely when selecting the required gas volume flow. Do not try to set gas volume flows between the numerical values. No gas will flow at interim settings.

"OFF"

"ON"

5.3.1.2 Continually adjustable flow outlet

- The flow gauge must be upright.
- Read the value on the upper edge of the float.

5.4 At the end of operation

- ➔ Manually close the valve on the compressed gas cylinder by turning it clockwise. Do not exert force.

Either use or completely empty the gas remaining in the pressure regulator (see Section "Safety information for working with gases"). The gauge must fall back to "0".

- ➔ The connection must be depressurised when the coupling is unfastened. It is strictly forbidden to unfasten the coupling when pressurised and using a tool except in an emergency. This will damage the gasket beyond repair, necessitating a subsequent replacement.

Downstream equipment and accessories can then be disconnected. Please also refer to the instructions for use issued by the manufacturer of the downstream equipment.

Use suitable protective caps to seal the connections on the compressed gas cylinder, the pressure regulator and any other appliances.

5.5 Changing the compressed gas cylinder

GB

- ➔ 1. Manually close the valve on the compressed gas cylinder by turning it clockwise. Do not exert force.
2. Either use or completely empty the gas remaining in the pressure regulator (see Section "Safety information").
The connection must be depressurised when the coupling is unfastened. It is strictly forbidden to unfasten the coupling when pressurised and using a tool except in an emergency. This will damage the gasket beyond repair, necessitating a subsequent replacement.
3. Downstream equipment can then be disconnected. Please also refer to the instructions for use issued by the manufacturer of the downstream equipment.
4. Threaded connections with hexagon nuts must be loosened using an appropriate tool. Hexagon nuts with right-handed threads (without a groove around them) must be loosened by turning them anticlockwise. Hexagon nuts with left-handed threads (with a groove around them) must be loosened by turning them clockwise (see Fig. 2).
5. Continue to Section "Connection to a compressed gas cylinder".

5.6 Decommissioning and storage

- ➔ 1. Manually close the valve on the compressed gas cylinder by turning it clockwise. Do not exert force.
2. Either use or completely empty the gas remaining in the pressure regulator (see Section "Safety information for working with gases"). The gauge must fall back to "0".
The connection must be depressurised when the coupling is unfastened. It is strictly forbidden to unfasten the coupling when pressurised and using a tool except in an emergency. This will damage the gasket beyond repair, necessitating a subsequent replacement.
3. Downstream equipment can then be disconnected. Please also refer to the instructions for use issued by the manufacturer of the downstream equipment.
4. Threaded connections with hexagon nuts must be loosened using an appropriate tool. Hexagon nuts with right-handed threads (without a groove around them) must be loosened by turning them anticlockwise. Hexagon nuts with left-handed threads (with a groove around them) must be loosened by turning them clockwise (see Fig. 2).
5. Use suitable protective caps to seal the connections on the compressed gas cylinder, the pressure regulator and any other appliances.
6. Store the pressure regulator in a clean, dust-free and dry place until the next time it needs to be used.
7. If the pressure regulator is out of use for a prolonged period, it should be cleaned (see Section "Maintenance and care").

Important!

It is also vital to ensure that maintenance intervals are observed for appliances put into storage. Only appliances which have been properly maintained whilst in storage may be put to use when taken out of storage.

6 MAINTENANCE AND CARE

6.1 Cleaning and disinfecting

Non-compliance with the maintenance and care instructions for the appliance can result in severe injury to the patient.

Foreign materials or fluids must not get inside the appliance. The pressure regulator may not be immersed in liquid!

The pressure regulator may not be sterilised!

6.1 Cleaning and disinfecting

- All of the inlet and outlet connections on the appliance must be sealed with suitable and clean protective caps/plugs prior to cleaning and disinfecting the equipment.

6.1.1 Cleaning

Wipe the surface of the pressure regulator with a dry disposable cloth or with a disposable cloth that has been moistened with clean water.

Cleaning agents containing alcohol or emollients may ignite if they come into contact with compressed oxygen or nitrous oxide.

6.1.2 Disinfection

Only use approved surface disinfectants. Observe the manufacturer's instructions for use.

Disinfectants that are based on the following materials are compatible:

- aldehydes
- quaternary ammonium compounds

Examples of approved preparations:

- Incidin rapid, Ecolab Deutschland GmbH
- terralin protect, Schülke & Mayr
- antifact extra, Schülke & Mayr

Disinfectants based on the following materials are not suitable as they may cause damage:

- halogen-splitting compounds
- strong organic acids
- oxygen-splitting compounds

GB

Users in the Federal Republic of Germany are recommended to always use disinfectants that are included in the latest VAH register (VAH: Verbund für Angewandte Hygiene e.V., / Association for Applied Hygiene, mhp Verlag GmbH, Wiesbaden). The VAH register also lists the active agents on which each disinfectant is based.

The recommendation of the active agents listed above applies in those countries that are unaware of the VAH register.

6.2 Performance test

A performance test must be carried out every time a cylinder is changed and at least every 6 months (if the appliance has been out of service).

Should faults be detected during performance testing or operation, the TESCOM EUROPE pressure regulator must be taken out of service immediately. Contact TESCOM EUROPE or an authorised specialist dealer to have the fault rectified.

A complete performance test includes:

- A visual inspection for mechanical damage
- Checking the system for leaks
- Checking the safety valve for leaks

6.2.1 Checking the system for leaks

1. Connect a filled compressed gas cylinder (see Section "Connection to a compressed gas cylinder").
2. Connect a suitable downstream appliance and start it up (see Section "Connection of downstream appliances").
3. Check that all screw connections and hose connections are secure.
4. Open the valve on the compressed gas cylinder slowly.
5. The cylinder pressure can now be read on the pressure regulator's inlet pressure gauge.
6. Close the cylinder valve again.
7. Watch the pointer on the inlet pressure gauge for approximately 1 minute. If the pointer remains constant, the system is leak-tight. If the pointer drops steadily, the system is not leak-tight. In this case you should proceed as follows:
 - 7.1 Mix a soap and water solution using perfumed soap.
 - 7.2 Open the compressed gas cylinder.
 - 7.3 Apply the solution to all screw connections. Bubbles will form wherever connections are not leak-tight.
 - 7.4 Close the compressed gas cylinder.
 - 7.5 Either use or completely empty the gas remaining in the pressure regulator (see Section "Safety information for working with gases"). The gauge must fall back to "0".

The connection must be depressurised when the coupling is unfastened. It is strictly forbidden to unfasten the coupling when pressurised and using a tool except in an emergency. This will damage the gasket beyond repair, necessitating a subsequent replacement.

7.6 In the event of leaks please replace damaged gaskets (see Section "Replacing gaskets") or contact TESCOM EUROPE.

7.7 Then check the system again to ensure that it is leak-tight.

If you are unable to eliminate the leak, the appliance must be repaired by TESCOM EUROPE or an authorised specialist dealer.

6.2.2 Checking the pressure relief valve for leaks

1. Connect a filled compressed gas cylinder (see Section "Connection to a compressed gas cylinder").
2. Connect a suitable downstream appliance and start it up (see Section "Connection of downstream appliances").
3. Check that all screw connections and hose connections are secure.
4. Open the valve on the compressed gas cylinder **slowly**.
5. Use your finger to check whether gas is escaping from the pressure relief valve.
6. If gas is escaping, you will need to take the pressure regulator out of service and have it repaired by the manufacturer (TESCOM EUROPE) or an authorised specialist dealer.

6.2.3 Replacing gaskets

1. It is absolutely essential that you wash your hands thoroughly before changing the gasket.
2. Take a small screwdriver or a similar tool and clean it thoroughly with washing-up liquid. Then rinse it with plenty of water and dry thoroughly. The tool must be clean and free from oil and grease in order to avoid the risk of explosion.

3. Using the clean tool, carefully remove the gasket.

Connection using the handwheel

- Make absolutely sure that the groove inside the gasket is not damaged when making the handwheel connection.
- If necessary, wipe the groove for the gasket clean with a dry cloth. You can also moisten the cloth with a little clean water. The safety and function of the appliance can only be assured if original TESCOM EUROPE spare parts are used.
- Carefully push the new gasket into the connection groove. Do not use a tool to do this under any circumstances.

Connection with pin index

- Make absolutely sure that the groove inside the gasket is not damaged.
- If necessary, wipe the sealing surface clean with a dry cloth. You can also moisten the cloth with a little clean water.
- The safety and function of the appliance can only be assured if original TESCOM EUROPE spare parts are used.
- Carefully push the new gasket onto the connection sealing surface. Do not use a tool to do this under any circumstances.

GB

Connection by means of a threaded connection


- Make absolutely sure that there is no damage to the groove or the connection piece sealing surface in which or on which the gasket is located.
- If necessary, wipe the groove for the gasket clean with a dry cloth. You can also moisten the cloth with a little clean water. The safety and function of the appliance can only be assured if original TESCOM EUROPE spare parts are used.
- Carefully push the new gasket into the connection groove. Flat sealing discs must be placed on the sealing surface of the connection piece for the pressure regulator. Do not use a tool to do this under any circumstances.

6.3 Maintenance and service life

The service life of the appliance is 12 years from the month and year of manufacture, which are stated on the type plate.

Fig. 5 Example of type plate

Year of manufacture Month of manufacture

The appliance must be disposed of once it has reached the end of its service life (see Section "Disposal").

Maintenance including - amongst other things - the replacement of all wear parts must be carried out after 6 years. Non-compliance with this deadline will invalidate the warranty provided by TESCOM EUROPE. The safety and function of the appliance can only be assured subject to maintenance being carried out by TESCOM EUROPE or an authorised company and provided that original TESCOM EUROPE spare parts are used.

Compressed gas cylinders:

Compressed gas cylinders must undergo repeat testing by a certified testing agency at regular intervals. Before using a compressed gas cylinder, check that the date for such repeat testing has not elapsed. Compressed gas cylinders whose date has elapsed must not be used. Repeat testing dates are stamped on compressed gas cylinders.

6.4 Warranty, customer service and dispatch

All products purchased from TESCOM EUROPE GMBH & CO. KG are covered by the statutory warranty from the date of purchase.

The warranty obligations and liability of TESCOM EUROPE GMBH & CO. KG for any consequences are excluded in the following cases:

- Operation after the end of service life
- Non-compliance with the instructions for use
- Use for any other than the intended purpose
- Operating errors
- Improper use or handling
- Negligence
- Tampering with the pressure regulator by unauthorised third parties, e.g. opening the appliance or replacing connection pieces
- Act of God
- Damage during transport caused by improper packaging
- Failure to maintain the appliance
- Failure to use original spare parts

Please also refer to the corresponding section in our terms of delivery on the delivery note.

You will find the manufacturer's address on the last page of these instructions for use.

Please contact our repairs department before returning any goods for repair or any faulty goods. In the interests of our staff we also require a declaration indicating whether – and if so which – harmful substances have been used with the appliance and which measures (e.g. rinsing) have been taken so that we can exclude any hazards during repair.

6.5 Troubleshooting

GB

Problem	Cause	Solution
Connection between cylinder valve and pressure regulator is not leak-tight.	Damaged gasket.	Replace gasket.
Outlet pressure rising, pressure relief valve blowing off.	Soiled or damaged valve seat.	Must be repaired by TESCOM EUROPE or an authorised specialist dealer.
Leaking around spring cap.	Faulty diaphragm.	Must be repaired by TESCOM EUROPE or an authorised specialist dealer.
Unusually high gas consumption.	System leaking.	Find and fix the leak. If the pressure regulator is leaking, it must be repaired by TESCOM EUROPE or an authorised specialist dealer.
Flow decreases over the duration.	Blocked inlet filter.	Must be repaired by TESCOM EUROPE or an authorised specialist dealer.
Flow outlet does not close.	Damaged flow outlet.	Must be repaired by TESCOM EUROPE or an authorised specialist dealer.
Flow outlet cannot be adjusted.	Damaged flow outlet.	Must be repaired by TESCOM EUROPE or an authorised specialist dealer.

6.6 Disposal

Please contact an approved and certified disposal company to ensure the proper disposal of the pressure regulator once it has reached the end of its service life. Your environmental protection representative or municipal authorities will be able to provide the relevant contact details. Pay attention to local, regional or domestic regulations and legislation.

7 COMBINATION WITH DOWNSTREAM EQUIPMENT

The user and operator are responsible for using and operating the pressure regulator and the equipment connected to it.

Downstream appliances and accessories must comply with the following criteria:

1. It must be a medical product with CE mark.
2. It must be marked with the same type of gas as the pressure regulator.
3. The downstream appliance must be suitable for the pressure regulator and the application. Details can be found in the technical specifications of the appropriate instructions for use.
Only use amagnetic accessories and appliances in the range of strong magnetic fields. Check that they are suitable prior to use.
4. Prior to use, the user must always check to ensure the material compatibility of the downstream equipment and appliances with the gas and/or gas mixture.

The temperature of the gas and the pressure regulator drops sharply during operation. This can produce hypothermia in the patient.

Should you have any questions or require specialist advice, please contact us or visit us on the Internet at www.tescom-europe.com.

You will find the address of the manufacturer of this pressure regulator on the last page of these instructions for use.

8 APPENDIX

8.1 Pressure regulator: Technical specifications

Subject to design changes.

8.1.1 Pressure regulator with pressure outlet

Permissible inlet pressure range: 1.1 to max. 30 MPa (11 to max. 300 bar) (depending on inlet connection)

Connection	Nominal outlet pressure *)	Nominal flow (air)
NFS 90-116	350±50 kPa (3.5±0.5 bar)	40 l/min
BS 5682	400±50 kPa (4.0±0.5 bar)	60 l/min
Other connections	450±50 kPa (4.5±0.5 bar)	80 l/min

*) at minimum inlet pressure (1.1 MPa) and nominal flow

Due to the design of the appliance, the outlet pressure increases linearly as the inlet pressure decreases (approx. 150 kPa (1.5 bar) at 30 MPa to 1.1 MPa (300 bar to 11 bar)).

The actual outlet pressure is determined by the flow rate and the instantaneous inlet pressure (see Section "Flow curves").

GB

The flow is also dependent on the type of gas. Please consult the table below for the nominal flow of gases other than air.

AIR	O ₂	N ₂ O	CO ₂	Xe	He
30 l/min	28 l/min	24 l/min	24 l/min	14 l/min	80 l/min
55 l/min	52 l/min	45 l/min	45 l/min	25 l/min	145 l/min
80 l/min	76 l/min	65 l/min	65 l/min	37 l/min	215 l/min

8.1.2 Pressure regulator flow outlet

Permissible inlet pressure range: 1.1 to max. 20 MPa (11 to max. 200 bar) (depending on inlet connection)

Configurable nominal gas volume flows and accuracy (gas temperature 20 ± 2 °C; 1013 hPa, free flow into the atmosphere)

Due to the design of the appliance, the gas volume flow varies in response to fluctuations in the pressure of the gas supply and temperature fluctuations, as well as when the flow resistance at the outlet changes, e.g. on account of downstream equipment.

8.1.2.1 Locking adjustable flow outlet

Flow setting range: 0.1 – 30 l/min (depending on the combination of orifices)

Accuracy when flow < 1.0 l/min: -20 / +15% of the displayed value

Accuracy when flow 1 - 30 l/min: ... ±15% of the displayed value

8.1.2.2 Continually adjustable flow outlet

Flow setting range: 0.027 – 32 l/min *) (depending on the gauge head)

Accuracy of appliances at max. flow ≤ 2 l/min: ±2.5% of the end scale value

Accuracy of appliances at max. flow > 2 l/min: ±10% of the displayed value

8.1.3 Test gas pressure regulator

Version	Inlet pressure range	Nominal outlet pressure *)	Nominal flow (air)
2.1 bar	0.52 to max. 20 MPa	210 ± 10kPa (2.1±1 bar)	30 l/min.
5 bar	1.1 to max. 20 MPa	500 ± 10kPa (5 ± 1 bar)	55 l/min.
7 bar	1.5 to max. 20 MPa	700 ± 10kPa (7 ± 1 bar)	100 l/min.

*) at minimum inlet pressure (0.52 / 1.1 / 1.5 MPa) and nominal flow

Due to the design of the appliance, the outlet pressure increases linearly as the inlet pressure decreases (approx. 150 kPa (1.5 bar) at 200 MPa to 1.1 MPa (20 bar to 11 bar)).

The actual outlet pressure is determined by the flow rate and the instantaneous inlet pressure (see Section "Flow curves").

The flow is also dependent on the type of gas. Please consult the table below for the nominal flow of gases other than air.

AIR	O ₂	N ₂ O	CO ₂	Xe	He
40 l/min	38 l/min	32 l/min	32 l/min	18 l/min	107 l/min
60 l/min	57 l/min	48 l/min	48 l/min	28 l/min	160 l/min
100 l/min	95 l/min	81 l/min	81 l/min	46 l/min	270 l/min

8.2 Pressure relief valve: Technical specifications

Opening pressure800±50 kPa (8±0.5 bar)
 Max. flow rate220 l/min. air

8.3 Permissible ambient conditions


Air humidity:5% to 95% rel. humidity, non-condensing

Air pressure

Operation:900 hPa to 1100 hPa
 Storage:115 hPa to 2000 hPa

Permissible temperature range for brass pressure regulators

Storage


Operation

Permissible temperature range for stainless steel pressure regulators

Storage

Operation


Weight:approx. 1.5 kg (depending on the model)


8.4 Flow curves (for AIR)

8.4.1 Pressure regulator with pressure outlet

GB

8.4.2 Test gas pressure regulator

8.5 Ordering information

See separate sheet D44698

9 SPECIAL INFORMATION

9.1 Information on the electronic high-pressure sensor


Compliance with relevant national safety regulations is mandatory when commissioning and operating this pressure transducer.

Non-compliance with the relevant regulations could result in severe injury and/or property damage. Only suitably qualified staff may work on these appliances.

Electrical connection



The pressure transducer is joined via connectors or cables to the electrical power supply. The correct terminal assignment can be seen on the following block diagrams. In addition, the terminal assignment, output signal and requisite auxiliary power are noted on the type plate of the pressure transducer.



An stable direct current voltage within the limits specified on the type plate is sufficient as an auxiliary power supply. In the case of pressure transducers with current output, care must be taken to ensure that the supply voltage is sufficient to guarantee the minimum requisite voltage for external display and analysis equipment.

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