Gebrauchsanweisung **D**
Instruction for Use **GB**
Instructions d'utilisation **F**

Druckminderer Wegamed Compact Wegamed Compact cylinder pressure regulator Régulateur-détendeur Wegamed compact

- NL** Deze gebruiksaanwijzing vindt u in uw landstaal op de bijgesloten cd-rom.
- E** Encontrará estas instrucciones de uso en su idioma nativo en el CD-ROM adjunto.
- I** Le presenti istruzioni per l'uso sono disponibili nella lingua del proprio Paese nel CD-ROM allegato.
- P** No CD incluído encontra estas instruções de uso na sua língua.
- S** Den här bruksanvisningen hittar du på ditt språk på den medföljande CD-ROM-skivan.
- N** Denne bruksanvisningen er tilgjengelig på ditt aktuelle språk på den vedlagte CD-ROM-en.
- FIN** Käyttöohje on suomenkielisenä laitteen mukana toimitettavalla CD-ROM-levyllä.
- DK** Denne brugsanvisning findes på vedlagte cd-rom på dit sprog.
- RUS** Данную Инструкцию по эксплуатации можно найти на своём родном языке на прилагаемом компакт-диске CD-ROM.

D INHALT

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1	For your safety and that of your patients	23	5	Maintenance and care ..	34
1.1	Safety information	23	5.1	Cleaning and disinfecting	34
1.2	Safety information for working with medical gases ..	24	5.2	Performance test	35
1.3	Rules and guidelines	25	5.2.1	Checking the system for leaks	35
1.4	Applicable standards	25	5.2.2	Checking the pressure relief valve for leaks	36
1.5	Symbols and information on the appliance.....	26	5.2.3	Changing gaskets	36
2	Product description	26	5.3	Maintenance	37
2.1	Intended purpose	26	5.4	Warranty, customer service and dispatch	37
2.2	Set-up and mode of operation.....	27	5.5	Troubleshooting	38
3	About the instructions for use.....	27	5.6	Disposal	38
3.1	Scope.....	27	6	Combination with downstream equipment	38
3.2	Target audience	27	7	Appendix	39
3.3	Improvements	27	7.1	Technical specifications.....	39
4	Installation and operation.....	28	7.2	Flow curves (air)	40
4.1	Transport and packaging	28	7.3	Classification.....	41
4.2	Connection to a compressed gas cylinder	28	7.4	UMDNS code.....	41
4.3	Connecting downstream equipment.....	32	7.5	Ordering information	41
4.4	Commissioning	33	7.6	Dimensions	42
4.5	At the end of operation.....	33			
4.6	Changing the compressed gas cylinder	34			
4.7	Decommissioning and storage	34			

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
1 FOR YOUR SAFETY AND THAT OF YOUR PATIENTS

1.1 Safety information

- Please read and observe the following information before using or performing maintenance on the TESCOM EUROPE WegaMed Compact.
- The WegaMed Compact may only be used for its intended purpose as described in these instructions for use (see Chapter 2.1 "Intended purpose").
- These instructions for use are integral to the cylinder pressure regulator and must be available to the user at all times.
- You must understand how the cylinder pressure regulator works before commencing operation with it.
- Never use an appliance if it is leaking or faulty!
- Operation of the WegaMed Compact enriches and where appropriate controls the corresponding gas in the patient environment. The WegaMed Compact may therefore only be used in accordance with a doctor's instructions.
- Medical gases are a medicinal product and as such may only be used in accordance with instructions issued by the doctor responsible for the care of the patient in question.
- Before commissioning the WegaMed Compact, it is absolutely vital to check that the gas type, the gas volume flow range and the connection are correct and that any upstream or downstream equipment has been fitted securely.
- Not suitable as a drive for surgical equipment.
- Not suitable for operation in the range of strong magnetic fields, e.g. magnetic resonance imaging.
- Only accessories bearing the CE mark that are compatible with the WegaMed Compact both in terms of intended use and connection may be used (see Chapter 6 "Combination with downstream equipment"). The use of third-party products can lead to functional failures and a lack of biocompatibility. Please note that in such cases any and all warranty and liability claims will be invalidated, if neither the accessories recommended in the instructions for use nor original spare parts have been used.
- To avoid misunderstandings, never leave the appliance unattended or connected to downstream equipment when it is not in use.
- Never try to open the appliance since this could result in the release of spring-loaded parts which might cause injury.

- To avoid the risk of infection or bacterial contamination, please read and observe the information in Chapter 5 "Maintenance and care".
- The WegaMed Compact must not be immersed in liquids or sterilised. Make sure that foreign bodies cannot get inside the appliance. This applies to dust and disinfectants in particular. Use the protective caps supplied.
- You must wash your hands before working on the oxygen supply.
- Ensure that the appliance cannot come into contact with oil or grease during operation or storage. Oil and grease as found for example in hand creams can react fiercely with gases under pressure (especially oxygen and nitrous oxide). Explosion hazard!
- Do not use antiseize or lubricants on the WegaMed Compact.
- Smoking and naked flames are strictly forbidden in the vicinity of valves and fittings carrying oxygen.
- Always proceed with caution when opening the cylinder valve.
- The WegaMed Compact connections must match those on the gas supply or downstream equipment. Adapters must not be used.
- The compressed gas cylinder must be secured at all times so that it cannot tip over.
- Make sure that the temperature of the gas administered corresponds to the patient's body temperature. Risk of injury to the patient.

1.2 Safety information for working with medical gases

Oxygen (O₂) increases combustibility and flammability as well as encouraging flammable materials to burn. Greases and oils (such as those found in hand creams or lubricants) can ignite on contact with oxygen under pressure. Good ventilation is absolutely essential. High levels of concentration can be damaging to health. Smoking and naked flames are forbidden.

Nitrous oxide (N₂O) has a strong anaesthetic effect. It increases combustibility and flammability as well as encouraging flammable materials to burn. Greases and oils (such as those found in hand creams or lubricants) can ignite on contact with nitrous oxide under pressure. High levels of concentration will introduce the risk of anoxia and apnoea. Good ventilation is absolutely essential. Smoking and naked flames are forbidden.

Carbon dioxide (CO₂) initially has an intoxicating effect; this is followed by a debilitating effect. Increased levels of concentration will introduce the risk of anoxia and apnoea. Good ventilation is absolutely essential.

Xenon (Xe) has a strong anaesthetic effect. High levels of concentration will introduce the risk of anoxia and apnoea. Good ventilation is absolutely essential.

High levels of **nitrogen (N₂)** can reduce O₂ content in breathable air in enclosed spaces. Good ventilation is absolutely essential when high levels of nitrogen are present.

1.3 Rules and guidelines

The following rules and guidelines must be observed from case to case, in addition to the aforementioned safety information:

1. Regulations for operators of medical equipment in Germany (Medizinprodukte-Betreiberverordnung or MPBetreibV for short)
2. Principles of accident prevention (BGV A1)
3. Material usage (BGR 500)
4. Information provided by the suppliers and manufacturers of compressed gases

1.4 Applicable standards

CE marking in accordance with EC Directive 2007/47/EC

The WegaMed Compact meets the requirements of the following standards:

- DIN EN ISO 10524
- CGA E-4
- UL 252

Optional inlet connection

- DIN 477: 2003
- DIN EN ISO 407: 2005
- BS341
- NF E 29-650
- NEN 3268
- ITC MIE AP7
- UNI 4406
- ISO 5145

Optional outlet connection

- EN 739 (NIST)
- DIN 13260
- NFS 90-116
- BS 5682
- UNI 9507
- SS8752430
- CGA V-5 (DISS)

1.5 Symbols and information on the appliance

The meaning of the symbols that are both attached to the appliance (e.g. type plate) and contained in the accompanying information, such as instructions for use or packaging, is explained in the following section.

GB

	<p>"Caution, consult the accompanying documentation!"</p> <p>This symbol indicates that particular attention must be paid to the relevant section or information to ensure the safe use of the product.</p>
	<p>"Manufacturer"</p> <p>The name of the manufacturer appears next to this symbol. It might be accompanied by the manufacturer's address.</p>
	<p>"Date of manufacture"</p> <p>The year of manufacture appears next to this symbol.</p>
	<p>"Free from oil and grease"</p> <p>This symbol denotes equipment that uses medical gases and must therefore be kept free from oil and grease.</p>
	<p>"Temperature restriction"</p> <p>This symbol indicates a temperature range, e.g. in connection with storage or operating conditions.</p>

2 PRODUCT DESCRIPTION

2.1 Intended purpose

The WegaMed Compact is used for a variety of medical equipment in the field of medicine and therapy to regulate the pressure at which medical gases are delivered. For this purpose, users connect it to compressed gas cylinders (up to 20 MPa = 200 bar cylinder pressure). It is not suitable for the operation of surgical equipment.

Reduction of 20 MPa (200 bar) inlet pressure maximum to a nominal outlet pressure of 450 kPa (4.5 bar).

Only to be used for the medical gas indicated on the WegaMed Compact:

- Oxygen (O₂),
- Nitrous oxide (N₂O),
- Carbon dioxide (CO₂),
- Air for artificial respiration (AIR),
- Xenon (Xe)

The regulator may only be used to supply suitable medical equipment in accordance with the technical specifications indicated on the type plate and in the instructions for use.

2.2 Set-up and mode of operation

The TESCOM EUROPE WegaMed Compact consists of a single-stage, direct-acting, spring-loaded pressure regulator with inlet pressure compensation. The TESCOM EUROPE WegaMed Compact is fitted with a pressure relief valve.

3 ABOUT THE INSTRUCTIONS FOR USE

3.1 Scope

These instructions for use are valid for the TESCOM EUROPE WegaMed Compact, part number D50770.

3.2 Target audience

These instructions for use are intended for anyone responsible for the installation, operation and monitoring of the above valves and fittings.

The operator has a duty to ensure that only personnel who are familiar with the general practice and safety measures associated with working with gases and gas valves and fittings, as well as with the fundamental regulations with regard to occupational safety and accident prevention, are permitted to work on the device. These personnel must also have been taught how to handle the pressurised device, and have read and understood the instructions for use.

It is absolutely essential that you read and understand at least Chapter 1 "For your safety and that of your patients" and Chapter 4 "Installation and operation" before attempting to operate the valves and fittings.

Operators must be given access to these instructions. Additional copies can be supplied on request.

3.3 Improvements

If you have any suggestions for improving the product or the instructions for use, we would be delighted to hear from you. Both our products and the instructions for use are subject to a process of continuous development and improvement. The address and telephone number for TESCOM EUROPE GMBH & CO. KG appear on the last page of these instructions for use.

4 INSTALLATION AND OPERATION

Responsibility for using and operating the pressure regulator and accessories lies with the user or operator.

GB

4.1 Transport and packaging

The TESCOM EUROPE WegaMed Compact may only be transported (e.g. delivery to the manufacturer) and stored in suitable, sturdy packaging (e.g. the original packaging). During transport and storage the inlet and outlet connectors must be covered with suitable protective caps.

4.2 Connection to a compressed gas cylinder

The use of dust-free, dry medical gases is essential if the TESCOM EUROPE WegaMed Compact is to function as intended.

Damaged or soiled appliances or gas cylinders must never be used. A check must be carried out to ensure that there is no evidence of damage or soiling on the connections on the gas cylinder and the WegaMed Compact. Damaged gaskets must be replaced. New gaskets can be purchased from TESCOM EUROPE or an authorised specialist dealer.

1. You must wash your hands thoroughly before commencing any work involving the oxygen supply. Hydrocarbon compounds (e.g. oils, greases, rubbing alcohols, hand creams or adhesive plaster) can cause explosive reactions if they come into contact with highly compressed gases.
2. Using the technical specifications on the type plate and in these instructions for use, check whether the pressure regulator in your possession is suitable for the intended purpose (type of gas, pressure, flow). The maximum permissible inlet pressure of the pressure regulator must be greater than or equal to the cylinder's filling pressure. Read the information provided by the supplier of the compressed gases carefully. Medical gases are a medicinal product and as such may only be used in accordance with instructions issued by the doctor responsible for the care of the patient in question.

3. The pressure regulator must undergo a performance test before every cylinder change. Should faults be detected during the performance test, the pressure regulator must be taken out of service. Contact TESCOM EUROPE or an authorised specialist dealer to have the fault rectified.

A complete performance test includes:

- A visual inspection for mechanical damage
- Checking the system for leaks (Chapter 5.2.1)
- Checking the pressure relief valve for leaks (Chapter 5.2.2)

4. Before connecting the TESCOM EUROPE WegaMed Compact, slowly open the valve on the compressed gas cylinder briefly in a well ventilated room or outdoors (this is to blow out any contamination).

Pressure surges caused by opening the valve too quickly can cause explosions! To avoid the risk of injury, the gas jet must be directed away from any persons in the vicinity.

GB

5. Remove the protective caps covering the connections on the pressure regulator and keep them in a safe place.
6. Screw the TESCOM EUROPE WegaMed Compact to the compressed gas cylinder. The connections must match - **do not use adapters!**

All connections must be clean and free from oil and grease. Do not use lubricants.

Handwheel connection

The pressure regulator can be fitted with a handwheel connection to make it easier to connect it to the cylinder valve. This connection is denoted by the knurled connection screw or nut and the elastomer gasket inside one of the grooves on the coupling connection.

When making the connection, it is absolutely essential to check the gasket for evidence of damage and soiling. If the gasket is damaged or soiled, the pressure regulator must not be used and has to be replaced. (Chapter 5.2.3 "Replacing gaskets"; ordering information appears in Chapter 7.5.)

Please note that this screw connection needs to be made without using tools.

As soon as the compressed gas cylinder is opened, the compressed gas pushes the gasket into the annulus between the valve on the compressed gas cylinder and the pressure regulator. This creates a reliable seal.

The connection must be depressurised when the coupling is unfastened. It is strictly forbidden to unfasten the coupling when pressurised and using a tool except in an emergency. This will damage the gasket beyond repair, necessitating a subsequent replacement.

Check that the pressure regulator and the compressed gas cylinder are securely connected before **slowly** opening the cylinder valve. If the connection is not leaktight, close the cylinder valve immediately.

Once the connection has been made correctly, run a performance test (see Chapter 5.2 "Performance test") to ensure safe and reliable operation.

Connection with PIN index

The pressure regulator can also be fitted with a PIN index. Please note that only the one original TESCO EUROPE gasket included in the scope of supply may be used.

GB

When making the connection, it is absolutely essential to check the gasket for evidence of damage and soiling. If the gasket is damaged or soiled, the pressure regulator must not be used and has to be replaced (Chapter 5.2.3 "Replacing gaskets"; ordering information appears in Chapter 7.5).

The pins sticking out of the PIN index must fit into the bore on the cylinder valve. The arrangement of these pins will vary depending upon the type of gas.

The entire PIN index - with the gasket - must slot into the sealing surface of the cylinder valve by hand (without having to exert excess force).

As soon as the compressed gas cylinder is opened, the compressed gas pushes the gasket into the annulus between the valve on the compressed gas cylinder and the pressure regulator. This creates a reliable seal.

Screw the tommy screw tight by hand (make sure that the tip of the tommy screw is sitting in the centre hole of the cylinder valve).

Check that the pressure regulator and the compressed gas cylinder are securely connected before **slowly** opening the cylinder valve. If the connection is not tight, close the cylinder valve immediately.

Once the connection has been made correctly, run a performance test (see Chapter 5.2 "Performance test") to ensure safe and reliable operation.

Responsibility for using and operating the pressure regulator and the equipment connected to it, as well as any accessories, lies with the user or operator.

4.3 Connecting downstream equipment

A variety of equipment can be connected downstream of the TESCO EUROPE WegaMed Compact.

Only suitable equipment may be connected to the pressure regulator. Chapter 6, "Combination with downstream equipment", describes permissible accessories.

When connecting equipment, please note the following:

1. Connections and gaskets must be clean and free from oil and grease.
2. There must be no evidence of damage on the gasket and sealing surfaces.

3. Correct and proper procedures must be followed when connecting the TESCOM EUROPE WegaMed Compact to downstream equipment. Please also read the instructions issued by the manufacturers of the downstream equipment you are using.
4. All downstream equipment, hoses and connecting elements must be marked with the same gas type.
5. Once the relevant connections have been made, a check must be carried out to ensure that all equipment is in working order.
6. Compliance with the Medizinprodukte-Betreiberverordnung (regulations for operators of medical equipment, Germany) is mandatory.

4.4 Commissioning

1. Connect and start up downstream equipment. Please also follow the instructions issued by the manufacturers of the downstream equipment you are using.

2. Open the valve on the compressed gas cylinder **slowly**. Open the handwheel on the cylinder valve as far as it will go and then turn it back half a revolution. The cylinder pressure will now be displayed on the inlet pressure gauge.

3. Read the cylinder pressure and determine the how long the gas in it will last. You might need to replace the compressed gas cylinder to ensure that the downstream equipment will remain in operation for the requisite time. Medical gases are a medicinal product and as such may only be used in accordance with instructions issued by the doctor responsible for the care of the patient in question.

The patient and the pressure indicated on the gauge must be monitored continuously.

If no gas is being drawn: close the valve on the compressed gas cylinder. The pressure regulator is not a check valve.

To avoid misunderstandings, never leave the appliance connected to the patient unsupervised or connected when it is not in use.

Make sure that gas from the patient is not able to get inside the pressure regulator (backflow). Should this happen, the appliance must be sent back to the manufacturer immediately for cleaning.

4.5 At the end of operation

Close the valve on the compressed gas cylinder by hand.

Use up the remaining gas in the pressure regulator or vent it fully. The gauge must fall back to "0". Wait until the pressure regulator has depressurised fully before unscrewing it - failure to do this can damage the gaskets.

Downstream equipment can then be disconnected. Please also refer to the instructions for use issued by the manufacturer of the downstream equipment.

Cover the connections with suitable protective caps.

4.6 Changing the compressed gas cylinder

1. Close the valve on the compressed gas cylinder by hand.
2. Use up the remaining gas in the pressure regulator or vent it fully. Compliance with the safety instructions for medical gases is mandatory. The gauge must fall back to "0". Wait until the pressure regulator has depressurised fully before unscrewing it - failure to do this can damage the gaskets.
3. Downstream equipment can then be disconnected. Please also refer to the instructions for use issued by the manufacturer of the downstream equipment.
4. Unscrew the pressure regulator from the valve on the compressed gas cylinder.
5. More information appears in Chapter 4.2 "Connection to a compressed gas cylinder".

4.7 Decommissioning and storage

1. Close the valve on the compressed gas cylinder by hand.
2. Use up the remaining gas in the pressure regulator or vent it fully. Compliance with the safety instructions for medical gases is mandatory. The gauge must fall back to "0". Wait until the pressure regulator has depressurised fully before unscrewing it - failure to do this can damage the gaskets.
3. Downstream equipment can then be disconnected. Please also refer to the instructions for use issued by the manufacturer of the downstream equipment.
4. Unscrew the pressure regulator from the valve on the compressed gas cylinder.
5. Cover the inlet and outlet connections on the TESCOM EUROPE WegaMed Compact with protective caps.
6. Store the pressure regulator in a clean, dust-free and dry place until the next time it needs to be used.
7. If the pressure regulator is out of use for a prolonged period, it should be cleaned (see Chapter 5 "Maintenance and care").

Important!

It is also vital to ensure that maintenance intervals are observed for appliances put into storage. Only appliances which have been properly maintained whilst in storage may be put to use when taken out of storage.

5 MAINTENANCE AND CARE

5.1 Cleaning and disinfecting

If exterior cleaning is required, please only use a clean cloth. The cloth can be dry or damp (with clean water). If, in exceptional cases, wet disinfection cannot be avoided, standard approved surface disinfectants must be used. Compliance with safety instructions for medical gases, as well as with guidelines for use, is mandatory.

Prior to cleaning or disinfection, all inlet and outlet connections on the TESCOM EUROPE WegaMed Compact must be fitted with suitable protective caps.

It is absolutely vital to ensure that neither cleaning agents nor disinfectants, nor other contaminants, are able to get inside the TESCOM EUROPE WegaMed Compact.

The appliance must not be immersed in liquids or sterilised.

5.2 Performance test

A performance test must be carried out every time a cylinder is changed and at least every 6 months (if the appliance has been out of service).

Should faults be detected during the performance test or when the appliance is in operation, the TESCOM EUROPE WegaMed Compact must be taken out of service immediately. Contact TESCOM EUROPE or an authorised specialist dealer to have the fault rectified.

A complete performance test includes:

- A visual inspection for mechanical damage
- Checking the system for leaks (Chapter 5.2.1)
- Checking the pressure relief valve for leaks (Chapter 5.2.2)

5.2.1 Checking the system for leaks

1. Connect a full compressed gas cylinder. The compressed gas connections must be screwed tight by hand only.
2. Connect a suitable item of downstream equipment and start it up.
3. Check that all screw connections and hose connections are secure.
4. Open the valve on the compressed gas cylinder slowly.
5. The cylinder pressure can now be read on the pressure regulator's inlet pressure gauge.
6. Close the cylinder valve again.
7. Watch the pointer on the inlet pressure gauge for approx. 1 minute. If the pointer remains constant, the system is leaktight. If the pointer drops steadily, the system is not leaktight. In this case you should proceed as follows:
 - 7.1 Mix a soap and water solution using perfumed soap.
 - 7.2 Open the compressed gas cylinder.
 - 7.3 Apply the solution to all screw connections. Bubbles will form wherever connections are not leaktight.
 - 7.4 Close the compressed gas cylinder.
 - 7.5 Use up the remaining gas in the pressure regulator or vent it fully. Compliance with the safety instructions for medical gases is mandatory. The gauge must fall back to "0". Wait until the pressure regulator has depressurised fully before unscrewing it - failure to do this can damage the gaskets.
 - 7.6 In the event of leaks please replace damaged gaskets (see 5.2.3 "Replacing gaskets") or contact TESCOM EUROPE.
 - 7.7 Then check the system again to ensure that it is leaktight.

If you are unable to make the system leaktight, repairs will be necessary.

5.2.2 Checking the pressure relief valve for leaks

GB

1. Connect a full compressed gas cylinder. The compressed gas connections must be screwed tight by hand only.
2. Connect a suitable item of downstream equipment and start it up.
3. Check that all screw connections and hose connections are secure.
4. Open the valve on the compressed gas cylinder slowly.
5. Use your finger to check whether gas is escaping from the pressure relief valve.
6. If gas is escaping, you will need to take the pressure regulator out of service and have it repaired by the manufacturer (TESCOM EUROPE) or an authorised specialist dealer.

5.2.3 Changing gaskets

1. It is absolutely essential that you wash your hands thoroughly before changing the gasket.

2. Take a small screwdriver or a similar tool and clean it thoroughly with washing-up liquid. Then rinse it with a generous amount of water and dry thoroughly.

The tool must be clean and free from oil and grease in order to avoid the risk of explosion.

3. Using the clean tool, carefully remove the gasket.

Handwheel connection

- **Important!** It is vital to take care to ensure that the groove inside the gasket is not damaged when making the handwheel connection.
- If necessary, wipe the groove for the gasket clean with a dry cloth. You can also moisten the cloth with a small amount of clean water. The safety and function of the appliance can only be assured if original TESCOM EUROPE spare parts are used.
- Carefully push the new gasket into the connection groove. Do not use a tool to do this under any circumstances.

Connection with PIN index

- **Important!** It is vital to take care to ensure that the groove inside the gasket is not damaged when making the PIN index connection.
- If necessary, wipe the sealing surface clean with a dry cloth. You can also moisten the cloth with a small amount of clean water.

The safety and function of the appliance can only be assured if original TESCOM EUROPE spare parts are used.

- Carefully push the new gasket onto the connection sealing surface.
Do not use a tool to do this under any circumstances.

5.3 Maintenance

Maintenance including - amongst other things - the replacement of all wear parts must be carried out every 6 years. Noncompliance with this deadline will invalidate the warranty provided by TESCOM EUROPE.

The safety and function of the appliance can only be assured subject to maintenance being carried out by TESCOM EUROPE or an authorised company and provided that original TESCOM EUROPE spare parts are used.

Compressed gas cylinders:

Compressed gas cylinders must undergo repeat testing (by the TÜV in Germany) at regular intervals. Before using a compressed gas cylinder, check that the date for such repeat testing has not elapsed. Compressed gas cylinders whose date has elapsed must not be used. Repeat testing dates are stamped on compressed gas cylinders.

5.4 Warranty, customer service and dispatch

- All products purchased from TESCOM EUROPE GMBH & CO. KG are covered by the statutory warranty from the date of purchase.
- TESCOM EUROPE GMBH & CO. KG will not honour this warranty in the event of (e.g.):
 - Noncompliance with the instructions for use
 - Maloperation
 - Improper use or handling
 - Negligence
 - Tampering with the appliance by unauthorised third parties, e.g. opening the appliance or replacing connection pieces
 - Act of God
 - Damage during transport caused by improper packaging
 - Failure to maintain the appliance
 - Failure to use original spare parts
- Please also refer to the corresponding section in our terms of delivery on the delivery note.
- You will find the manufacturer's address on the last page of these instructions for use.
- Please contact our repairs department before returning any goods for repair or any faulty goods. In the interests of our staff we also require a declaration indicating

whether – and if so which – harmful substances have been used with the appliance and which measures (e.g. rinsing) have been taken to exclude any hazards during repair.

GB

5.5 Troubleshooting

Problem	Cause	Solution
Connection between cylinder valve and pressure regulator is not leaktight.	Damaged gasket.	Replace gasket.
Outlet pressure rising, pressure relief valve blowing off.	Soiled or damaged valve seat.	Send to manufacturer or authorised specialist dealer for repairs.
Leaking around spring cap.	Faulty diaphragm.	Send to manufacturer or authorised specialist dealer for repairs.
Unusually high gas consumption.	System leaking.	Find and fix the leak. If the Wegamed Compact is not leaktight it should be sent to TESCOM EUROPE GMBH & CO. KG or an authorised specialist dealer for repairs.

5.6 Disposal

Please contact an approved and certified disposal company to ensure the proper disposal of the Wegamed Compact. Your environmental protection representative or municipal authorities will be able to provide the relevant contact details.

6 COMBINATION WITH DOWNSTREAM EQUIPMENT

- Responsibility for using and operating the pressure regulator and the equipment connected to it lies with the user or operator.
- Please refer to the technical specifications in these instructions for use for information about the outlet pressure and flow rate of the pressure regulator; compliance with these specifications is mandatory.

Should you have any questions or require specialist advice, please contact us or visit us on the Internet at www.tescom-europe.com

7 APPENDIX

7.1 Technical specifications

Permissible inlet pressure range (p3 to p1)

Air (AIR)	1.1 to 20 MPa (11 to 200 bar)
Oxygen (O ₂)	1.1 to 20 MPa (11 to 200 bar)
Xenon (Xe)	1.1 to 6 MPa (11 to 60 bar)
Nitrous oxide (N ₂ O)	1.1 to 6 MPa (11 to 60 bar)
Carbon dioxide (CO ₂)	1.1 to 6 MPa (11 to 60 bar)

Nominal outlet pressure (p2)

set to a fixed value of 1.1 MPa (11 bar) inlet pressure p1

Connection	Nominal outlet pressure p2	Nominal flow rate
NFS 90-116	350±50 kPa (3.5±0.5 bar)	40 l/min.
BS 5682	400±50 kPa (4.0±0.5 bar)	40 l/min.
Other connections	450±50 kPa (4.5±0.5 bar)	40 l/min.

The nature of the appliance design is such that the outlet pressure increases as the inlet pressure (cylinder pressure) decreases (approx. 50 kPa (0.5 bar) at 200 MPa to 1.1 MPa (200 bar to 11 bar)). The actual outlet pressure is determined by the flow rate and the instantaneous inlet pressure.

Pressure relief valve


Opening pressure (8±0.5 bar)
Max. flow 220 l/min. air

Operating environment


Air humidity: 5% to 95% rel. humidity, non-condensing

Temperature

Storage:

Operation:


GB

Weight: 800 g (varies according to model)

Subject to design changes.

7.2 Flow curves (air)

The curve shown relates to operation with an outlet connection. Running both outlet connections at the same time significantly reduces the outlet pressure at each outlet connection.

7.3 Classification

Pursuant to EC Directive 2007/47/EC
Appendix IX Class IIb

7.4 UMDNS code

Nomenclature for medical equipment 13-323 high-pressure gas regulators


7.5 Ordering information

Connection configuration	Type of gas	Replacement gasket for cylinder connection
DIN 477 Part 1	O ₂	D50810
DIN 477 Part 1 (cylinders as of 3 l)	N ₂ O	D41414
DIN 477 Part 1 (cylinders up to 3 l)	N ₂ O	D43921
DIN 477 Part 1	AIR	D50809
DIN 477 Part 1	CO ₂	D50809
DIN 477 Part 1	He	D50809
DIN 477 Part 1	Xe	D50809
DIN EN ISO 407 (PIN index)	O ₂	D45574
DIN EN ISO 407 (PIN index)	N ₂ O	D45574
DIN EN ISO 407 (PIN index)	AIR	D45574
DIN EN ISO 407 (PIN index)	CO ₂	D45574
DIN EN ISO 407 (PIN index)	He	D45574
DIN EN ISO 407 (PIN index)	Xe	D45574
BS 341	AIR	1129023
BS 341	O ₂	1129023
BS 341	N ₂ O	D44195
BS 341	Xe	1129023
NF E 29-650	AIR	D44212
NF E 29-650	O ₂	D44212
NF E 29-650	N ₂ O	D44212
NEN 3268	AIR	D50809
NEN 3268	O ₂	D41414
NEN 3268	N ₂ O	D50809
ITC MIE	O ₂	D44212
UNI 4406	O ₂	D50809

7.6 Dimensions

Dimensions can vary from model to model. All dimensions are approximate measurements in millimetres.

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