

Branson™ Sonifier Accessories

Support Stand with Rod

Overview

The support stand with stainless steel rod accommodates Branson's converter and accessories. The converter clamp supports the ultrasonic stack (converter and horn) and is easily adjusted to properly position the horn in the sample.

The stand features a black enameled cast-iron base with zinc-plated rod. A full acoustic sound enclosure is also available for sound-sensitive environments.



Stand Measurements

Model 8800	
Stand Base	9" l x 6"w x 24" h (229 x 140 mm)
Rod	24"h x 1/2" dia. (609mm x 13mm)
Branson Part No. 101-063-1108	Stand with rod
Branson Part No. 101-063-1109	Converter clamp for use with Branson SFX150 Sonifiers
Branson Part No. 101-063-1110	Converter clamp for use with Branson SFX250, SFX550 Sonifiers

For more information:
www.Emerson.com/Branson

BRANSON™



EMERSON™