

Product Origin List

Today's Date: March 9, 2016

Date Last Updated: March 1, 2016

This matrix is designed to communicate the HTS, ECCN and *possible* origin of MMI equipment. Due to the complex nature of our manufacturing facilities, it is not possible to provide the exact origin of multiple-origin products without serial number information. See Foot Note 1 for specific origin information.

Country Code Descriptions

CN = China

DE = Germany

DK = Denmark

GB = Great Britain

MX = Mexico

NL = Netherlands

RO = Romania

SG = Singapore

US = United States

Model	Product Type	Origin	HTS	ECCN
Sensors				
CMFS-Series	CMFS007M, CMFS025M, CMFS040M, CMFS050M, CMFS075M, CMFS100M, CMFS150M	US	9026.90.2000	EAR99
	CMFS010M, CMFS015M	US or RO	9026.90.2000	EAR99
	CMFS025K, CMFS040K, CMFS050K, CMFS075K, CMFS100K, CMFS150K, CMFS010H, CMFS015H, CMFS025H, CMFS050H, CMFS100H, CMFS150H, CMFS010P, CMFS015P, CMFS025P, CMFS050P, CMFS100P, CMFS150P	US	9026.90.2000	EAR99
CMF-Series	CMF010, CMF025, CMF050, CMF100	MX	9026.90.2000	EAR99
	CMF200	MX	9026.90.2000	EAR99
	CMF200M	CN or MX	9026.90.2000	EAR99
	CMF300	MX	9026.90.2000	EAR99
	CMF300M	CN or MX	9026.90.2000	EAR99
	CMF350M, CMF350A	MX	9026.90.2000	EAR99
	CMF400	MX	9026.90.2000	EAR99
	CMF400M	CN or MX	9026.90.2000	EAR99
	CMFHC2A, Y CMFHC2G, M	US US or CN	9026.90.2000	EAR99
	CMFHC3A CMFHC3G, M CMFHC3Y	US US or CN US	9026.90.2000	EAR99
	CMFHC4M	US or CN	9026.90.2000	EAR99
D-Series	DH100, DH150, DH300	US	9026.90.2000	EAR99
	DS300, DS600	US	9026.90.2000	EAR99
	DL065, DL100, DL200	US	9026.90.2000	EAR99
	DT065, DT100, DT150	US	9026.90.2000	EAR99
CNG-Series	CNG050	CN, MX, RO	9026.90.2000	EAR99

Category	Product Type	Origin	HTS	ECCN
F-Series	F025H, F025A, F025B, F025P F050H, F050A, F050B, F050P F100H, F100A, F100B F200H	MX	9026.90.2000	EAR99
	F100P	MX	9026.90.2000	EAR99
	F025S, F050S, F100S, F200S	MX, CN, RO	9026.90.2000	EAR99
	F300	MX or CN	9026.90.2000	EAR99
H-Series	H025, H050, H100, H200, H300	MX	9026.90.2000	EAR99
Low Flow-Series	LF2M, LF3M, LF4M	US	9026.90.2000	EAR99
R-Series	R025S, R050S, R100S, R200S	MX, CN, RO	9026.90.2000	EAR99
T-Series	T025, T050, T075, T100, T150	US	9026.90.2000	EAR99
Density /Viscosity-Series	7826, 7828	MX	9026.80.2000	EAR99
	7827, 7829	MX	9027.80.4590	EAR99
	3098, 7812, 7835, 7845, 7847	GB	9026.80.2000	EAR99
	782781, 782782, 782784	GB	9027.90.5495	EAR99
	782791	MX	9027.90.5495	EAR99
	782992, 782993, 782994	GB	9026.90.5495	EAR99
	SGM, GDM	RO	9026.80.2000	EAR99
	CDM	US	9026.80.2000	EAR99
	FDM	MX	9026.80.2000	EAR99
	FVM, HFVM	MX	9027.80.4590	EAR99
Transmitters and Peripherals				
	1500, 2500	MX	9026.90.2000	EAR99
	1700, 2700	US, MX, CN, NL, RO	9026.90.2000	EAR99
	5700	US, CN, NL	9026.90.2000	EAR99
	2200S	MX	9026.90.2000	EAR99
	2400S	MX or CN	9026.90.2000	EAR99
	3100, 3300, 3350, 3500, 3700	MX, CN	9026.90.2000	EAR99
	7950, 7951, 7965	GB	9026.10.2040	EAR99
	FC475	DE	9031.80.8085	5A992
	775 Remote Mount	SG	8517.62.0050	5A991
	775 Integrally Mounted onto Remote Transmitter	Origin of Transmitter	9026.90.2000	EAR99
	775 Integrally Mounted onto Transmitter/Sensor	Origin of Sensor	9026.80.2000	EAR99
	1420	US	8517.62.0050	5A991
	9739MVD	MX	9026.90.2000	EAR99
	LFT1, LFT3,LFT4, LFT6	US	9026.90.2000	EAR99

Category	Product Type	Origin	HTS	ECCN
	LFT2, LFT5, LFT8	MX	9026.90.2000	EAR99
	PLK	US	9026.90.2000	EAR99
	FMT	MX	9026.90.2000	EAR99
	RFT9739	MX	9026.90.2000	EAR99
	BARRIERKITCSAE, BARRIERKITATEXE, BARRIERKITULE	DK	9032.89.2000	EAR99
	PUCK (Core Processor)	MX or CN	9026.90.2000	EAR99
Cable				
	CFEP, CLFT, CMVD, CPLT	DE	8544.11.0050	EAR99
Rosemount Magnetic Flowmeters & Vortex Flometers				
Go to	http://rmintra.emersonprocess.com/sites/rmt/importexport/ProductAttributeInformation/Rosemount_Product_Matrix.pdf			

Note 1: For products with multiple origins, the origin is determined by the following serial number ranges:

Sensors:

US Origin: 10000000 - 10999999 and 12000000 - 12999999
 Mexican Origin: 11000000 - 11999999 and 14000000 - 14999999
 Chinese Origin: 13000000 - 13999999
 GB Origin: 260000 - 299999
 Romanian Origin: 15000000 - 15999999

Transmitters:

Mexican Origin: 9000000 - 9999999 and 25000000 - 25999999 and 0392300000 - 0392999999 and 2182800000 - 2182999999 and 2294800000 - 2294999999
 Thailand Origin: 16000000 - 16999999 and 27000000 - 34999999
 Chinese Origin: 33000000-34999999 and 13000000 - 13999999
 US Origin: 10000000 - 10999999 and 12000000 - 12999999
 Romanian Origin: 15000000 - 15999999
 Netherlands Origin: 19000000 - 19999999

Note 2: The HTS number for a flowmeter (sensor and transmitter), where the sensor and transmitter are intended for use together, (whether or not physically attached) is 9026.80.2000. The Country of Origin for a flowmeter (sensor and transmitter), where the sensor and transmitter are intended for use together, is that of the sensor.

Note 3: For RMA and Restock sensors, see the above serial number ranges for country of origin.