

PRODUCT NOTIFICATION

ROSEMOUNT GAS CHROMATOGRAPHS

July 19, 2016

**TO: GC Global Sales Channel**

**SUBJECT: Product Life Cycle Status of the Model 1000 Gas Chromatograph**

The Model 1000 gas chromatograph became “OBSOLETE”, per the Rosemount Analytical Product Life Cycle Policy, July 2013 after 5 years of being “INACTIVE”. As the product has been in “OBSOLETE” status for over 2 years, availability of spares and service is only on a “best effort” basis. The recommendation is to replace the Model 1000 when repairs become necessary.

The successors to the Model 1000 are the 700XA and 1500XA. The 700XA and 1500XA incorporate many new features, including the new XA series of analytical valves, an easy to maintain electronics section, standard redundant Ethernet connections, an optional Foundation fieldbus interface, historical chromatogram storage and interoperability with MON 2020.

Rosemount Analytical recommends customers purchase a new 700XA for applications that require certification and metrology approvals and where space may be limited. It recommends customers purchase a new 1500XA for applications that may require an air bath oven, or the additional column and valve capacities of a traditional process gas chromatograph oven. Additional information on the 700XA can be found at [700XA product information website](#). Additional information on the 1500XA can be found at: [1500XA product information website](#).

Regards,

Bonnie Crossland  
Product Marketing Manager  
Gas Chromatographs  
ROSEMOUNT