

American Recovery and Reinvestment Act Compliance Certification

Rosemount certifies that the following instrumentation products are manufactured in the United States, in compliance with §1605 of the American Recovery & Reinvestment Act, as amended by Interim Rule Subpart 25.6, Volume 74 Fed. Reg. 14623-33 (Mar. 31, 2009) and Final Rule FAC 2005-45, Volume 75, No 167 (Effective Date Oct. 1, 2010):

- Rosemount Model 8800 Vortex Flowmeter
- Rosemount Model 3051 Pressure Transmitter
- Rosemount Model 2051 Pressure Transmitter
- Rosemount Model 3051S (“A1020” option code) Pressure Transmitter
- Rosemount Model 3051SMV Multivariable Transmitter
- Rosemount Model 3051SFA ProBar Flowmeter
- Rosemount Model 3051SFC Compact Orifice Flowmeter
- Rosemount Model 3051SFP Proplate Flowmeter
- Rosemount Model 3095 Multivariable Transmitter
- Rosemount Model 3095MFA Mass Probar Flowmeter
- Rosemount Model 3095MFC Compact Orifice Mass Flowmeter
- Rosemount Model 3095MFP Mass Proplate Flowmeter
- Rosemount Model 4600 Pressure Transmitter
- Rosemount Model 405 Compact Orifice
- Rosemount Model 485 Annubar
- Rosemount Model 1195 Integral Orifice
- Rosemount Model 285 Annubar
- Rosemount Model 1595 Orifice Plate
- Rosemount Model 1199 Remote Seal
- Rosemount Model 3300 (R3616 option code) Guided Wave Radar Level Transmitter
- Rosemount Model 5300 (R3616 option code) High Performance Guided Wave Radar Level Transmitter
- Rosemount Model 0078 Platinum Wire-Wound Temperature Sensor (Excluding High Temperature Options)
- Rosemount Model 0068Q Series 68Q Sanitary RTD Temperature Sensor
- Rosemount Model 0058C Series 58C Cut-to-Fit RTD Sensor
- Rosemount Model 0183 Thermocouple Temperature Sensor



Matthew LaJoy
Contracts Negotiator

Rev. E (November 1, 2010)

1 November, 2010

Date