

REGISTRATION APPROVAL

Wednesday, December 22, 2010
Emerson Process Management, Rosemou
2400 Barranca Parkway
Irvine, CA
92606

ATTENTION : Dana Crowley

Our File : 8340 [0 F

With reference to your submission respecting the registration of the item below, for legal use in the province, please note we have surveyed, approved and registered this design as noted.

MANUFACTURER : Emerson Process Management, Rosemount Analytical	
CATALOG OR DRAWING : 23311-00 Rev. R	
ITEM : Toroidal Conductivity Sensor	CRN : OF12064.53

We wish to point out that every fitting must be constructed strictly in accordance with the registered design.

Fitting registrations are required to be resubmitted for validation after ten (10) years from the registration date in accordance with CSA B51, Clause 4.2.1. The date of expiry for this registration is February 11, 2018.

Sincerely,



Shane Kleisinger

Codes and Standards Compliance

REMARKS: CRN registered under reciprocal agreement. Design code: Manufacturer's standards supported by proof test - ISA 82-03 Table 18-1. Scope of Registration: Model 228-02,-03-220PSI & Model 228-04,-05-150PSI.