

ROSEMOUNT GAS CHROMATOGRAPHS

July 19, 2016

TO: GC Global Sales Channel**SUBJECT: Product Life Cycle Status of the Model 1000A Gas Chromatograph**

The Model 1000A gas chromatograph is being moved to “OBSOLETE” status, per the Rosemount Analytical Product Life Cycle Policy, effective August 1, 2016. Emerson will make efforts to support the gas chromatograph for the next 2 years. However, due to the age and the speed of obsolescence of electronic components, it is difficult to obtain all of the necessary spare parts. Some parts are already obsolete and when current stocks run out no more parts will be available.

The successors to the Model 1000A are the 700XA and 1500XA. The 700XA and 1500XA incorporate many new features, including the new XA series of analytical valves, an easy to maintain electronics section, standard redundant Ethernet connections, an optional Foundation fieldbus interface, historical chromatogram storage and interoperability with MON 2020.

Rosemount Analytical recommends customers purchase a new 700XA for applications that require certification and metrology approvals and where space may be limited. It recommends customers purchase a new 1500XA for applications that may require an air bath oven, or the additional column and valve capacities of a traditional process gas chromatograph oven. Additional information on the 700XA can be found at [700XA product information website](#). Additional information on the 1500XA can be found at: [1500XA product information website](#).

Regards,

Bonnie Crossland
Product Marketing Manager
Gas Chromatographs
ROSEMOUNT