

ROSEMOUNT™ 3051S MULTIVARIABLE™ TRANSMITTER

Measure differential pressure, static pressure and process temperature using a single device.

Simplify Mass Flow Measurement

To achieve the most accurate flow measurement, differential pressure, static pressure, and process temperature must be accounted for.


- Measuring only differential pressure can result in >6% error due to density changes.
- Adding static pressure measurement helps reduce error to <3% but temperature still affects mass flow.
- The addition of process temperature can reduce error down to <1%.

Taking each one of these three measurements separately involves complex calculations and multiple devices.

Accurate Mass Flow Calculations

The Rosemount™ 3051S MultiVariable™ Transmitter measures differential pressure, static pressure and process temperature in a single device while dynamically calculating mass flow 22 times per second within the transmitter.

The transmitter's ability to internally calculate fully compensated mass flow improves accuracy and reduces measurement uncertainty and complexity. This multivariable design also reduces pipe penetrations, impulse piping, and connection systems for additional cost savings.



Rosemount 3051S MultiVariable Transmitter

Key Features

- Delivers mass flow at .65% flow accuracy
- Dynamically corrects for all variables
- Compatible with all fluids, primaries, AGA, ISO and ASME Standards

For more information, visit Emerson.com
or contact your local Emerson Sales Representative



ROSEMOUNT 3051S MULTIVARIABLE TRANSMITTER


Fully Compensated Mass Flow Within the Transmitter

Rosemount 3051SMV Measures:

Differential Pressure
Static Pressure
Process Temperature

Enter:

Process Fluid
Primary Element
Line Size



Rosemount 3051SMV Calculates:


Density
Velocity
Viscosity
Beta Ratio
Gas Expansion
Discharge Coefficient
Velocity Approach
Reynolds Number

Receive:

Mass Flow
Energy Flow
Totalized Flow
Static Pressure
Differential Pressure
Process Temperature
Volumetric Flow
Compressibility Factor

Reduce Measurement Uncertainty

The Rosemount 3051SMV compensates for pressure and temperature variation, greatly reducing the measurement uncertainty associated with DP only measurement.



Steam: Normal Flow, 13.4 lb/s, 8% pressure change, 5% temperature change, Cost \$10/1,000 lbs
Natural Gas: Normal Flow, 1.8 lb/s, 14% pressure change, 10% temperature change, Cost \$5/1000 ft3
Air: Normal Flow, 1.3 lb/s, 13% pressure change, 11% temperature change, Cost \$.20/1000 ft3

For more information, visit
[Emerson.com/Rosemount3051SMV](https://www.emerson.com/Rosemount3051SMV)

The Emerson logo is a trademark and service mark of Emerson Electric Co.
 Rosemount is a registered trademark of one of the Emerson family of companies.
 All other marks are the property of their respective owners.
 © 2025 Emerson Electric Co. All rights reserved.

00807-0400-4803 Rev AB